



Multi Drug Rapid Test 7 Parameter



Specification

- Application** : For the simultaneous detection of amphetamine, methamphetamine, tetrahydrocannabinol, cocaine, barbiturates benzodiazepine & morphine/opiates in human urine.
- Sensitivity** : ≥ 1000 ng/ml of Ampethamine (AMP)
 ≥ 1000 ng/ml of Methamphetamine (MET)
 ≥ 50 ng/ml of Tetrahydrocannabinol (THC)
 ≥ 300 ng/ml of Cocaine (COC)
 ≥ 300 ng/ml of Benzodiazepine (BZO)
 ≥ 300 ng/ml of Morphine/Opiates (MOP)
 ≥ 2000 ng/ml of Carisoprodol (CAR)
- Storage & Stability** : The sealed pouches in the test kit and the kit components may be stored between 2 - 30°C till the duration of the shelf life as indicated on the pouch/carton. DO NOT FREEZE.

 catalogue



AKL 20103420324