




Multi Drug Rapid Test 8 Parameter



Specification

Application : For the simultaneous detection of amphetamine, methamphetamine, tetrahydrocannabinol, cocaine, barbiturates benzodiazepine, morphine/opiates & synthetic cannabinoids in human urine.

Sensitivity : ≥ 1000 ng/ml of Amphetamine (AMP)
 ≥ 1000 ng/ml of Methamphetamine (MET)
 ≥ 50 ng/ml of Tetrahydrocannabinol (THC)
 ≥ 300 ng/ml of Cocaine (COC)
 ≥ 300 ng/ml of Benzodiazepine (BZO)
 ≥ 300 ng/ml of Morphine/Opiates (MOP)
 ≥ 2000 ng/ml of Carisoprodol (CAR)
 ≥ 50 ng/ml of Synthetic cannabinoids (K2)

Storage & Stability : The sealed pouches in the test kit and the kit components may be stored between 2 - 30°C till the duration of the shelf life as indicated on the pouch/carton. DO NOT FREEZE.



AKL 20103420350